

Product sheet – Chrysal Professional 2 Premium T-Bag

Chrysal Professional 2 Premium T-Bag is a conditioning product for all cut flowers (with the exception of Anthurium), intended for use by wholesalers, bouquet makers and retailers.

It contains ingredients for the hydration of flowers and includes the correct amount of food needed during transport & retail phases.

Chrysal Professional 2 Premium T-Bag is part of the Chrysal guarantee concept.



General

- Conditioner for all cut flowers, with the exception of Anthurium.
- Suitable for use in cold storage, transport and retail phase.
- Contains the correct amount of nutrients for the natural development of the flower, without stimulating full development during this phase.
- Available in T-Bag formats for 1 L, 2 L and 4 L of solution, loose or continuous string.

Effects

- Keeps flowers beautiful up to 50% longer compared to water alone.
- Stimulates water uptake, thus preventing “bent-neck”.
- Lowers the pH of the water.
- Keeps flowers, especially roses, and foliage in optimum condition.
- Delays first stem death of roses.
- Better colour retention of the flowers.
- Leaves and stems remain firm and green, no slimy stems.
- Improves vase life of flowers compared to water alone.
- Contents of T-Bag dissolve clear and without odour in water.

Applications

- Suitable for all types of tap water.
- Can be dosed by hand or automatically.
- Visual dosage control evident.

Savings

- Reduces flower wastage.
- Reduces water consumption and labour costs, as the solution can be re-used.
- Extends distribution and sales period.
- Reduces the number of possible consumer claims.



Recommendations

- Store under dry, cool and dark conditions, preferably between 5°C and 25°C.
- Shelf life: 18 months in sealed packaging and stored under the right conditions.
- PH of the ready to use solution: 3.5 – 5.5 depending on the alkalinity of the water.
- Dosage: 1 T-Bag per volume of water as indicated on the T-Bag.
- Solution can be re-used for a maximum of 5 days, depending on temperature and use.
- Do not mix residual solutions with freshly made ones.
- Do not use in unprotected metal buckets or vases (zinc, copper, iron).
- Residual solutions may be drained into the sewer system.